

Bi 1x

E. coli Growth Curves

1 Overview

In this lab, you will investigate growth of the bacterium *E. coli*, watching the growth in two ways. First, you will determine growth rate in a solution of media by measuring optical density over time. The optical density is proportional to the mass of bacteria in the solution, within a certain dynamic range. You will measure the growth rate of *E. coli* using different sugars to compare growth rates. You will also repeat and expand upon a foundational experiment in molecular biology by growing *E. coli* in media with different sugars and observing the characteristics of growth curves.

Second, you will watch individual bacteria grow and divide using fluorescence microscopy. In this case, the mass of bacteria is proportional to the fraction of your images occupied by bacteria.

Through these two different approaches, you will learn about the rates that bacterial populations convert media into new life. You will also learn about the ability of bacteria to respond to changes in external conditions by turning on a suite of diverse catabolic enzymes.

2 Background

2.1 Bacterial growth

As you have already seen in our introductory experiments, bacteria can grow and divide very rapidly. A solution containing approximately 100 bacteria can grow to well over 10^9 bacteria—roughly equal to the entire human population—overnight. As with human or animal populations growing under optimal conditions, the rate of population increase is proportional to the current population. This means that

$$\frac{dN}{dt} = rN, \tag{1}$$

$$N(t) = N_0 e^{rt}, \tag{2}$$

where $N(t)$ is the number of bacteria in your sample at time t and $N_0 = N(0)$. The growth rate r depends on many features of the environment, such as the growth medium, temperature, and presence of antibiotics. We will investigate how the sugar present in the media changes characteristics of bacterial growth.

2.2 The Beer-Lambert Law and optical density

To determine the density (concentration) of bacteria in liquid media, we load a plate reader with the suspension of cells. The plate reader is a parallel spectrophotometer, in which light with a wavelength of 600 nm enters the sample at an initial intensity (I_0). Cells, having a different index of refraction as the surrounding medium, randomly reflect and scatter light out of the incident light path. The amount of scattering is proportional to the mass of cells at any given point in the sample. (Caveat: to get an accurate reading, cell density must be low enough that light is only scattered once in its journey through the sample. Often, this means that cells must be diluted to get an accurate reading.) Upon exiting the sample, the light intensity is reduced by the cells from I_0 to $I = I_0 e^{-\epsilon lc}$, where ϵ is the wavelength-dependent molar absorption coefficient of the species in question, l is the path length (that is, the total length of the trajectory of a light beam through the sample), and c is the concentration. This is known as the Beer-Lambert Law. The transmission efficiency of the sample is the ratio of the final and initial light intensity, $T = I/I_0$, and the reported quantity is called the absorbance,

$$A = -\ln T = \ln \frac{I_0}{I} = \epsilon lc. \quad (3)$$

Note that in order to compute the absorbance, I_0 must be measured. This is done by measuring a **blank** sample, which contains only the media in which the bacteria are grown, but no bacteria. Internally, the instrument has a value for I_0 , which we will call I_0^{inst} . So, the absorbance value given by the instrument for the measured absorbance is

$$A^{\text{inst}} = \ln \frac{I_0^{\text{inst}}}{I}. \quad (4)$$

Measurement of the blank gives

$$A_{\text{blank}}^{\text{inst}} = \ln \frac{I_0^{\text{inst}}}{I_0}, \quad (5)$$

since $I_0 = I_{\text{blank}}$ because there are no bacteria to scatter light. Therefore, the absorbance is

$$A = A^{\text{inst}} - A_{\text{blank}}^{\text{inst}} = \ln \frac{I_0^{\text{inst}}}{I} - \ln \frac{I_0^{\text{inst}}}{I_0} = \ln \frac{I_0}{I}. \quad (6)$$

The general rule of thumb is that the accurate absorbance readings lie in the range of $0.01 < A < 1$. (Serial dilutions of your sample may be necessary to bring the absorbance into this range.) In the context of bacterial growth, the absorbance A is often called **optical density**, abbreviated OD_{600} for absorbance of light at a wavelength of 600 nm.

2.3 Calibration of optical density

In the previous section, we established that the optical density is proportional to cell concentration. To determine the constant of proportionality (the product ϵl), we dilute a sample of known OD to a very low concentration of bacteria. We then spread a small volume of this diluted sample onto a plate with a good growth medium. We incubate the plate and count the number of bacterial colonies present. Provided the initial solution was dilute enough, each colony should correspond to a single bacterium. Since we know how much we diluted the original sample, we can then back-calculate its concentration. This gives us the constant or proportionality between OD_{600} and cell concentration. A good rule of thumb is that an OD_{600} of 0.1 corresponds to about 10^8 cells/mL.

2.4 Your Mission

- Determine bulk growth rates and diauxic transitions using different carbon sources.
- Determine the correspondence between OD_{600} and cell concentration in liquid medium.

3 Protocol

3.1 Bulk Assay

Prior to your involvement, overnight cultures of wild type *E. coli* will be inoculated in both LB medium and M9 minimal medium with glucose in a shaker at 37°C. These cultures will be completely saturated and at stationary phase.

3.1.1 Bulk Assay OD Calibration

You will measure OD_{600} using the lab plate reader.

1. Pipette 1 mL of overnight culture into an Eppendorf tube.
2. Prepare a series of eight ten-fold dilutions of your culture (i.e., 10^0 to 10^{-7}) in 500 μ L volumes using the provided Eppendorf tubes. Note that because you may transfer some of one dilution to the next, your final volume may not be exactly 500 μ L for some of the dilutions. For example, say you have made 500 μ L of a 10^{-3} dilution and you choose to make your 10^{-4} dilution by mixing 50 μ L of your 10^{-3} dilution with 450 μ L of media, then you would end up with 450 μ L of your 10^{-3} dilution.
3. Add 300 μ L of each dilution to the wells assigned to you in the 96 well plate.
4. Once all groups have loaded their samples, measure the OD in the plate reader and write down the values you measured for your samples.

5. Using the multichannel pipette, prepare a drip plate using the dilutions from the well plate.
6. To do this, use a multichannel pipette to remove 10 μL from each of your dilutions from the 96-well plate. Dispense onto the top of agar plate, fully expelling the volume.
7. Gently tilt the plate until vertical, allowing the droplets to streak to the opposite edge of the plate. Avoid mixing the droplets.
8. Place the plate back flat onto the bench top until the liquid dries.

This will be done one group at a time, so while you are waiting, prepare the media as described in the next section.

3.1.2 Bulk Assay Growth Preparation

1. Each group will be working with glucose and one additional sugar, chosen from xylose, sorbitol, galactose, rhamnose, ribose, and arabinose. You will be given 50 μL of a sugar at 20% concentration, 10 μL of glucose at 20% concentration, and 5 mL of M9 minimal media. In one 1.5 mL Eppendorf tube, make 1 mL of M9 minimal media with 0.2% glucose. In a 15 mL Falcon tube, make 5 mL of M9 minimal media at 0.2% of your other provided carbon source.
2. In three 1.5 mL Eppendorf tubes, prepare 1 mL of the sugar:glucose solutions as worked out in the prelab. Make sure to label your tubes.
3. Pipette 300 μL of each of your 4 solutions (sugar alone and 3 sugar:glucose mixes) into 3 wells of your assigned row (so you have 3 replicates for each). Note down the order in your lab notebook.
4. Using a multichannel pipette, pipette 3 μL of the 1:1000 M9 minimal media culture dilution that your TAs have prepared from a clean well boat into your assigned row. Make sure to fully expel the cells.

After all groups have filled their respective row, you will discuss and decide together how to fill the remaining row. Here you can get creative! Would you like to mix two different carbon sources other than glucose? Mix at different ratios? Use lower or higher sugar concentrations?

Your TAs will prepare 3 replicate wells of *E. coli* in glucose and M9 minimal media the same way in which you prepared your assigned sugar. They will also prepare 3 replicate wells of *E. coli* in LB (a media rich with nutrients for bacteria). You will use these results to compare to the growth characteristics in your assigned sugar.

Also note, as mentioned above, that having a blank is essential for OD₆₀₀ calculation. Your TAs will also include 3 minimal media blanks and 3 LB blanks in the top row of the 96-well plate.

The 96-well plate will then be placed in the plate read, kept at a constant temperature of 37°C, and OD₆₀₀ measurements will be made for the next 24 hours as bacteria grows.




Figure 1: Layout of samples in the 96-well plate. Ratios in row B through G are all sugar:glucose. In row H, you are encouraged to try any sugar combination that you fancy.

4 Assignment

Problem 0 (Summary).

Write a summary of this experiment, its objectives, and conclusions between one paragraph and one page in length.

Problem 1 (Serial dilution ([prelab](#))).

Write a sensible protocol for the serial dilutions for the OD calibration.

Problem 2 (Bulk growth recipes ([prelab](#))).

You will make two M9 minimal media solutions, one with 0.2% glucose and another with 0.2% of an assigned sugar. Fill out the table below with volumes to bring solutions to desired ratios.

Table 1: Media recipes.

sugar:glucose	0.2% sugar in M9	0.2% glucose in M9	Total volume
4:1			300 μ L
2:1			300 μ L
1:1			300 μ L

Problem 3 (OD calibration).

Count the colonies on your plates and use this value to determine a calibration factor between OD_{600} and absolute number of bacteria. You should also assess for what range of concentrations is the OD linear in the number of cells. How does this calibration factor compare to our rule of thumb that $OD_{600} = 0.1$ corresponds to 10^8 cells/mL?

Problem 4 (Plots of bulk growth curves: All the sugars).

In one plot, plot all growth curves obtained from pure sugar sources. That is, plot all curves obtained from your session's plate in wells containing glucose + minimal media, xylose + minimal media, galactose + minimal media, etc. Additionally, plot curves obtained from LB. Comment on what you see.

Problem 5 (Plots of bulk growth curves: Glucose and assigned sugar).

- For each of the three glucose + minimal media wells, plot OD_{600} versus time. All three curves should be on the same plot. (You already did this in the previous problem, but make a fresh plot.)
- Perform a regression to obtain growth rates for each of the three replicates. Report the growth rates and add the best fit lines to the plot.
- Repeat parts (a) and (b) the three wells containing your sugar + minimal media and also for the LB wells. All nine curves should be on the same plot for ease of comparison.
- How do the growth rates compare? Why do you think this is?

Problem 6 (Bulk growth: measuring growth rates before and after diauxic shifts).

For your three ratios (glucose:assigned sugar as 1:1, 1:2, and 1:4), plot your OD_{600} measurements versus time. First, note what you see qualitatively: where does the shift in growth rates occur for each ratio, and why do you think that is?

Next, for each ratio, select regions of the growth curve before and after the diauxic shift. Perform a regression to get the growth rates separately for each region.

How do these growth rates compare with the ones you fit in the previous problem for *E. coli* grown in glucose and your assigned sugar alone? Is this what you expected? Explain what may be happening in the cells to allow for the observed shift in growth rates.

Problem 7 (Data and code).

Attach all code used for data analysis. Be sure to write down the exact serial dilution protocols you did to calibrate your OD₆₀₀. This will likely be what you did in prelab question 1, but make sure you note what you actually did in lab.